





ECG Simulator

ST-6



Main features

- Fast and easy to use
- Low cost
- Heart beat signal
- Low battery signal
- Low power
- More than 100 hours autonomy (using alkaline battery)
- Light 170 grams (without battery)





Technical Features

Heart Beat	
Range	40 bpm aprox. - 120 bpm
Precision	+/- 2% (para 120 bpm)
Optical indicators	
Heart beat (QRS) / Discharged battery	1 Led bicolor GREEN/RED
Keys	
	Switch ON, switch OFF. Increased of the Amplitude signal x5
Safety	Cl. IP
In Accordance with the standard	72/23/EEC; EN61010 (IEC1010)
Power Supply	
Type	1 alkaline battery 6LR61 or rechargeable battery
Power consumption	6 mA
Autonomy	> 100 h. with an alkaline battery 30 h. with the rechargeable battery (175 mAh)
Keyboard	Membrane
Weight	170 grams (without battery)
Dimensions	142 x 85 x 40 mm

Standard accessories:

- 1 Alkaline battery 9V 6LR61
- 1 ST-6 User Manual

Optional Accessories:

- Carriage case
- Rechargeable battery 9Volt
- Battery Charger 9Volt

International:
E-mail: export@stelec.com
Web: www.stelec.com

Tel.: +34 93 5862747
Fax: +34 93 5862749

ST-Electromedicina, s.a.
Atenas, 23-25 Nave 2
Pol. Ind. Cova Solera
08191-RUBÍ
(Barcelona) SPAIN